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
解析モジュール	パラメータ	Acquisition/Replay での演算方法	Review での演算方法
BP	+dP/dt	Calculated between two successive end of peak detection times	Calculated between the end diastolic point and the systolic point
	-dP/dt	Calculated between the point that clears Minimum Pulse Height to the end of peak detection time	Calculated between the peak and the end of the peak detection time
LVP	Min	Calculated between the previous cycles Min Slope and the current cycles LVEDP point	Calculated between the Min Slope mark and the following cycle's LVEDP mark
ECG	PH	P direction is determined by the shape of the P wave relative to the T-P segment	P direction is based on the shape of the P wave relative to the points at which the P start and P end are marked
PAF	TV	Reports the maximum volume over the entire cycle	Reports the volume at the start of expiration as the Tidal Volume
PCR	TV	Reports the maximum volume over the entire cycle	Reports the volume at the start of expiration as the Tidal Volume
PT	Avg	Calculated between the point after clearing Peak Validation time for the previous cycle to the point at which peak validation time is cleared for the current cycle	Calculated from the current cycle's Minimum mark to the point prior to the following cycles Minimum mark
	Rate	Calculated based on the interval between the completion of the previous cycles peak validation time and the current cycles peak validation time	Calculated off the duration between the previous cycle's Maximum mark and the current cycles Maximum Mark
	+d_/dt	Calculated between the Start Contraction Mark and the completion of Peak Validation Time	Calculated between the Minimum Mark and the Maximum Mark

	-d_/dt	Calculated between the Maximum Mark and the point that clears the minimum pulse height	Calculated between the Maximum Mark and the next cycle's Minimum Mark
	+dd/dt	Calculated between the Start Contraction point and the peak validation time	Calculated between the Minimum Mark and the Maximum Mark
	-dd/dt	Calculated between the maximum point and the point that clears the minimum pulse height	Calculated between the Maximum Mark and the next cycle's Minimum Mark
URP	MV	Averaged value calculated off of the averaged TV and averaged BPM values when running in a multiple epoch logging rate or time-based logging rate	Averages are calculated from individual epochs, cycles, being averaged
	MVadj	Averaged value calculated off of the averaged TVadj and averaged BPM values when running in a multiple epoch logging rate or second logging rate	Averages are calculated from individual epochs, cycles, being averaged

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